



PACEY CUFF™
URETHRAL CONTROL DEVICE

Pacey Cuff is a urethral control device, designed for men who want to take control of their urinary continence care.

DR. JACK PACEY, MD FRCSC

URINARY INCONTINENCE

25 MILLION
AMERICANS
SUFFER

TREATABLE IN
MAJORITY OF
PEOPLE

AVAILABLE
SOLUTIONS
WERE OUTDATED

So we created a device that works painlessly and lets you get back to your active lifestyle by stopping urinary leakage without restricting blood flow

FEATURES

WHY THE PACEY CUFF?

1 BLOCK URINE
URETHRAL PRESSURE PAD BLOCKS UNWANTED URINE FLOW

2 CIRCULATION
WINDOWED HOOD ENSURES NORMAL BLOOD FLOW

3 COMFORTABLE
VELCRO STRAP AND COMPRESSION TAB HELPS ACHIEVE BEST FIT

4 QUALITY
PROFESSIONALLY DESIGN ALLOWS 6-8 MONTH USE

	Pacey Cuff	Disposable Pads	Other Clamps
	Proven to significantly reduce or stop leaks caused by Stress Urinary Incontinence (SUI) while maintaining blood flow in the penis	Absorbs leaks rather than stopping them	Stop leaks but also restrict blood flow to the penis, potentially causing pain
	Can use all day* (simply loosen to void, then re-tighten). Recommend loosening hourly for a brief period.	Must be changed when wet / saturated	MUST be released every hour to avoid compression injury / circulation compromise
	Reusable - lasts 6 to 8 months	Disposable - thrown away every few hours	Reusable
	Compact and discreet	Bulky and messy	Can be bulky or heavy
	Comfortable when fitted properly	Uncomfortable and awkward	Must be released hourly , can be painful
	Better for the environment	Single Use: Ends up in landfill daily - 500 years to decompose	Better for the environment
	Less than £150 a year (if used daily)	£200-£500	-

GET YOUR CUFF TODAY




website: MyPelvicHealth.co.uk
phone: +44 1923 23 77 95
email: info@mypelvichealth.co.uk

STEP 1

Circumference

Measurement is done mid-shaft around a non-erect (flaccid) penis.



Penis Circumference


SIZE	INCHES	CENTIMETERS
S	1 - 2.75"	3 - 7 cm
M	2.75 - 4	7 - 10 cm
L	4 - 5.5"	10 - 14 cm

* if your measurement is between two sizes we recommend choosing the larger size.

STEP 2

Minimum Penis Length

Measurement is from base to tip of non-erect penis.



Minimum Penis Length

SIZE	INCHES	CENTIMETERS
S	1.25"	3 cm
M	1.75"	4 cm
L	2"	4.5 cm