



afex Active – Contents:	Variant – Brief Size (above Hip)	NHS Code:
1 x High Receptacle	Small (28-30")	A100-S-08
1 x Briefs (Variable Size)	Medium (32-34")	A100-M-08
1 x Re-usable 500 ml Collection Bag	Large (36-38")	A100-L-08
1 x Cleansing Container	XL (40-42")	A100-XL-08
1 x Cleansing Solution – 118ml	2X (44-46")	A100-2X-08
1 x Instructional DVD	3X (48-50")	A100-3X-08

afex Core Supporter – Contents:	Variant – Brief Size (above Hip)	NHS Code:
1 x High Receptacle	Small (28-30")	A100S-5105-C
1 x Core Briefs (Variable Size)	Medium (32-34")	A100M-5105-C
1 x Re-usable 500 ml Collection Bag	Large (36-38")	A100L-5105-C
1 x Cleansing Container	XL (40-42")	A100XL-5105-C
1 x Cleansing Solution – 118ml	2X (44-46")	A1002X-5105-C
	3X (48-50")	A1003X-5105-C

afex Sedentary – Contents:	Variant – Brief Size (above Hip)	NHS Code:
1 x Low Receptacle	Small (28-30")	A110-S-08
1 x Briefs (Variable Size)	Medium (32-34")	A110-M-08
1 x Re-usable 1500 ml Collection Bag	Large (36-38")	A110-L-08
1 x Cleansing Container	XL (40-42")	A110-XL-08
1 x Cleansing Solution – 118ml	2X (44-46")	A110-2X-08
1 x Instructional DVD	3X (48-50")	A110-3X-08

afex Mobility Assisted – Contents:	Variant – Brief Size (above Hip)	NHS Code:
1 x Low Supersoft Receptacle	Small (28-30")	A120-S-08
1 x Briefs (Variable Size)	Medium (32-34")	A120-M-08
1 x Re-usable 1200 ml Collection Bag	Large (36-38")	A120-L-08
1 x Cleansing Container	XL (40-42")	A120-XL-08
1 x Cleansing Solution – 118ml	2X (44-46")	A120-2X-08
1 x Tube Assembly Complete	3X (48-50")	A120-3X-08

afex Urine Collection Bags (30-40 day lifespan with continuous use)	NHS Code:
500 ml Bag	A400-B-08
500 ml Vented Bag	A400-V-08
1200 ml Leg Bag	A400-L-08
1500ml Extra Capacity Bag	A400-E-08
2000ml Bed Bag	A400-X-08

afex Receptacles (3-4 month lifespan with continuous use)	NHS Code:
High Receptacle	A300-F-08
Low Receptacle	A300-L-08
Low Supersoft Receptacle	A300-S-08

afex Extension Tubes and Stabilizer Systems	NHS Code:
Tube Assembly Complete	A700-A-08
Tube Stabilizer	A700-S-07
Collection Bag Stabilizer (for Extra Capacity Bag)	A710-S-09

afex Leg Bag Holder (Lifespan 3-4 months)	NHS Code:
Afex Leg Bag holder (Black)	A800-B-07
Afex Leg Bag holder (White)	A800-W-07

afex Urine Collection (thigh) Bags (30-40 day lifespan with continuous use)	NHS Code:
250ml Bag	A400S
500 ml Bag	A400-B-08
500 ml Vented Bag	A400-V-08
1500ml Extra Capacity Bag	A400-E-08

afex Briefs – Size (Above Hip)	afex Briefs (White)	afex Briefs (Black)	afex Core Supporter Briefs (White)
Extra Small (24-26")	A200-XS-08	A200-XS-08	A200-CS-XS
Small (28-30")	A200-S-08	A200-S-08	A200-CS-S
Medium (32-34")	A200-M-08	A200-M-08	A200-CS-M
Large (36-38")	A200-L-08	A200-L-08	A200-CS-L
XL (40-42")	A200-XL-08	A200-XL-08	A200-CS-XL
2X (44-46")	A200-2X-08	A200-2X-08	A200-CS-2X
3X (48-50")	A200-3X-08	A200-3X-08	A200-CS-3X

**NEW**  
High Pubic Pressure  
Afex Core



### iMEDicare afex HOME Visit Service

If you identify a patient who you think will benefit from an afex day or night system, please call iMEDicare on 01923 23 7795. We will be delighted to conduct a home assessment and fit a customised trial kit, providing comprehensive patient review and feedback. afex can now cater for all male continence scenarios 24/7 from mild stress incontinence to the full force urination associated with urgency, in a very cost-effective manner.

UK and ROI Distributor:  
iMEDicare Ltd, Unit 11, Shakespeare Ind. Est.  
Shakespeare Street, Watford WD24 5RR

CE EC REP: Eurolink (Europe) Ltd, Avalon House  
Marcham Road Abingdon, Oxon OX14 1UD, UK  
Manufacturer: Arcus Medical, 2327 Distribution Street  
Charlotte, NC 28203 USA



is the answer when  
Nature Calls..

On NHS FP10  
Prescription

Comfortable, Convenient and Cost-effective

## Fights Odour, Rash and Embarrassment



**Comfortable:**  
Loose fit, adhesive free receptacle with double-wall design reduces exposure to residual urine.



**Discreet:**  
Briefs conceal system under clothing. Double-sealed collection bag effectively contains odor.



**Convenient:**  
Collection bag connects simply with our twist-n-lock design and easily drains with our one-touch port.

Up to 50% cheaper than Pads



+44(0)1923 23 77 95  
contact@imedicare.co.uk



iMEDicare  
Pelvic Health Naturally...



**Adhesive-Free EXternal**  
Collection System for Men

CE  
Manufactured in the USA  
by Arcus Medical LLC

[www.mypelvichealth.co.uk](http://www.mypelvichealth.co.uk)



# 1 āfex<sup>®</sup> Active



The āfex Active system can be used for any general out-door activity. Unretracted Penis context. It has three main components: an ergonomic High receptacle, boxer briefs (size determined by hip measurement) and a 500 ml collection bag. Possible Bag Sleeve accessory not included.



# 2 āfex<sup>®</sup> Sedentary



The āfex Sedentary is supplied with an extra-capacity bag (1500 ml), variable waist-size boxer brief (96% cotton/ 4% Lycra) low style receptacle (Super-soft option available) and is better suited to long-term sitting positions.



# 3 āfex<sup>®</sup> Mobility Assisted



Featuring a low supersoft receptacle for comfortable extended daily wear that accommodates a wide range of motion, āfex Mobility Assisted, enables wheel-chair bound and mobility assisted men to lead active out-door lifestyles.



# 4 āfex<sup>®</sup> Night Guard



For Nocturnal applications, all āfex Night Guard can be adapted for those who sleep in a fixed, stable position by using a high or low style receptacle (depending on body shape and penile length) connected to a 2-litre bed bag via an extension tube. For lower volume applications, the extra capacity bag can be positioned between the legs.



# 5 āfex<sup>®</sup> Core Supporter



āfex Core Supporter is ideal for either high activity or retracted Penis contexts. It has three main components: an ergonomic High receptacle, Core Supporter jock strap style briefs that provide high pubic pressure (size determined by hip measurement) and a 500 ml collection bag. Possible Bag Sleeve accessory (not included).



The āfex system makes managing chronic male incontinence easier than ever. The 96% cotton boxer briefs discreetly and securely hold a plastic (PVC) curved receptacle (3-4 months lifespan) inside a front pocket that channels urine into a 500 ml reusable collection bag (12 ml vinyl - 4 to 6 week lifespan). Uncomfortable tubes or straps are not required - just attach the bag directly to the receptacle through the briefs for hours of use without maintenance or changing, at very low cost.

## Why choose āfex?

**Pubic Pressure** to hold receptacle against the body wall comes from a tight but comfortable stretch āfex Brief (4% Lycra)

**Secure Positioning** of the penis centrally within the receptacle by padded secure wrap hugging the shaft of the penis at the open end of the receptacle. Wrap absorbs splash-back and prevents "popping out" with active movement




**Double Walled** inner receptacle lining drains urine away from the penis, helping to keep skin dry



Comfortable 96% Cotton Boxer Briefs

Shape Retaining Latex-Free Receptacle

Direct-Connect Vinyl Bag



## āfex<sup>®</sup> User Success

requires 1) Flaccid Penis to dangle over the testicles when seated 2) Dexterity 3) Cognitive Skills



- ✓ No Adhesives Used
- ✓ Receptacle is Breathable with Natural Relaxed Fit
- ✓ Complete System is Latex Free
- ✓ Simple Three Part System Comes Ready to Use
- ✓ Double-Wall Receptacle Keeps Skin Away From Urine
- ✓ Receptacle with Direct-Connect Bag is Fast Flow Capable
- ✓ Pressure-Free Collection Bag with Activity Tabs

## When to use āfex?

All Types of Male Urinary Incontinence - Stress - Urge - Mixed - Overflow - Functional incontinence - Pad saturating in public - Pad odour issue - Adhesive from penile sheath causing soreness / excoriation