

AFEX Starter Kit Prescription Advice

Hospital:.....

Dear Doctor:..... Date:..... Signed.....

Regarding patient:.....

Name Printed.....

TITLE.....

Department.....

ADDRESS.....

.....

.....

Postcode.....

Surname:	Hospital No.:
First names(s):	DOB:
Address:	
SPACE FOR PATIENT IDENTIFICATION LABEL	

Your patient was seen in our clinic today. Following physical assessment and full instruction in correct usage it has been established that the AFEX starter kit and accessories recommended below will prove most suitable and effective for him as an incontinence management solution. Can you please prescribe:

afex Active – Contents:	Variant – Brief Size (above Hip)	NHS Code:	Quantity
1 x High Receptacle	Small (28-30")	A100-S-08	
1 x Briefs (Variable Size)	Medium (32-34")	A100-M-08	
1 x Re-useable 500 ml Collection Bag	Large (36-38")	A100-L-08	
1 x Cleansing Container	XL (40-42")	A100-XL-08	
1 x Cleansing Solution – 118ml	2X (44-46")	A100-2X-08	
1 x Instructional DVD	3X (48-50")	A100-3X-08	

AFEX INCONTINENCE MANAGEMENT SYSTEMS - Starter Kits

Please tick or circle starter kit required as determined by patients above hip size - secure fit recommended to enhance receptacle pubic pressure

afex Core Supporter – Contents:	Variant – Brief Size (above Hip)	NHS Code:	Quantity
1 x High Receptacle	Small (28-30")	A100S-5105-C	
1 x Core Briefs (Variable Size)	Medium (32-34")	A100M-5105-C	
1 x Re-useable 500 ml Collection Bag	Large (36-38")	A100L-5105-C	
1 x Cleansing Container	XL (40-42")	A100XL-5105-C	
1 x Cleansing Solution – 118ml	2X (44-46")	A1002X-5105-C	
	3X (48-50")	A1003X-5105-C	

afex Sedentary – Contents:	Variant – Brief Size (above Hip)	NHS Code:	Quantity
1 x Low Receptacle	Small (28-30")	A110-S-08	
1 x Briefs (Variable Size)	Medium (32-34")	A110-M-08	
1 x Re-useable 1500 ml Collection Bag	Large (36-38")	A110-L-08	
1 x Cleansing Container	XL (40-42")	A110-XL-08	
1 x Cleansing Solution – 118ml	2X (44-46")	A110-2X-08	
1 x Instructional DVD	3X (48-50")	A110-3X-08	

DELIVERY:

Patients are requested to take their prescription to their local community pharmacy who should order directly on

Tel: 01923 23 7795
or Fax:01923 804 206
orders@imedicare.co.uk

or via 'Specials' depts of wholesalers.

FREE delivery to the Pharmacy from iMEDicare LTD.

iMEDicare LTD,
Unit 11 Shakespeare Industrial Estate,
Shakespeare Street, Watford WD24 5RR
www.iMEDicare.co.uk
Tel 01923 23 7795
contact@imedicare.co.uk

afex Mobility Assisted –Contents:	Variant – Brief Size (above Hip)	NHS Code:	Quantity
1 x Low Super-soft Receptacle	Small (28-30")	A120-S-08	
1 x Briefs (Variable Size)	Medium (32-34")	A120-M-08	
1 x Re-useable 1200 ml Collection Bag	Large (36-38")	A120-L-08	
1 x Cleansing Container	XL (40-42")	A120-XL-08	
1 x Cleansing Solution – 118ml	2X (44-46")	A120-2X-08	
1 x Tube Assembly Complete	3X (48-50")	A120-3X-08	

AFEX Accessories Prescription Advice

afex Briefs – Size (Above Hip)	afex Briefs (White)	Quantity	afex Briefs (Black)	Quantity	afex Core Supporter Briefs (White)	Quantity
Extra Small (24-26")	A200-XS-08		A200-XS-08		A200-CS-XS	
Small (28-30")	A200-S-08		A200-S-08		A200-CS-S	
Medium (32-34")	A200-M-08		A200-M-08		A200-CS-M	
Large (36-38")	A200-L-08		A200-L-08		A200-CS-L	
XL (40-42")	A200-XL-08		A200-XL-08		A200-CS-XL	
2X (44-46")	A200-2X-08		A200-2X-08		A200-CS-2X	
3X (48-50")	A200-3X-08		A200-3X-08		A200-CS-3X	

afex Urine Collection Bags (30-40 day lifespan with continuous use)	NHS Code:	Quantity
250 ml Bag	A400S	
500 ml Bag	A400-B-08	
500 ml Vented Bag	A400-V-08	
1200 ml Leg Bag	A400-L-08	
1500ml Extra Capacity Bag	A400-E-08	
2000ml Bed Bag	A400-X-08	

afex Receptacles (3-4 month lifespan with continuous use)	NHS Code:	Quantity
High Receptacle	A300-F-08	
Low Receptacle	A300-L-08	
Low Supersoft Receptacle	A300-S-08	

afex Extension Tubes and Stabilizer Systems	NHS Code:	Quantity
Tube Assembly Complete	A700-A-08	
Tube Stabilizer	A700-S-07	
Collection Bag Stabilizer (for Extra Capacity Bag)	A710-S-09	

afex Leg Bag Holder (Lifespan 3-4 months)	NHS Code:	Quantity
Afex Leg Bag holder (Black)	A800-B-07	
Afex Leg Bag holder (White)	A800-W-07	

AFEX INCONTINENCE MANAGEMENT SYSTEMS - Starter Kits

Please tick or circle quantity of accessories required (For briefs - measure patients above hip size - secure fit recommended).

DELIVERY:

Patients are requested to take their prescription to their local community pharmacy who should order directly on

Tel: 01923 23 7795
or Fax: 01923 804 206
orders@imedicare.co.uk

or via 'Specials' depts of wholesalers.

FREE delivery to the Pharmacy from iMEDicare LTD.

iMEDicare LTD,
Unit 11 Shakespeare Industrial Estate,
Shakespeare Street, Watford WD24 5RR
www.iMEDicare.co.uk
Tel 01923 23 7795
contact@imedicare.co.uk